



Be Right™



GANIMEDE Automatic Laboratory Analyser

Analysis unit with integrated digestion for serial laboratory analysis .

The GANIMEDE system consists of three components: the sample changer, the control unit and the analysis units. The analysis unit with integrated digestion is optionally available for the determination of ortho and total phosphate or for the determination of total nitrogen (TN_b). The compact system makes it possible to work flexibly at different measuring stations. The cordless control unit communicates with one or both analysis units. Sample data is acquired via the barcode scanner or a PC and processed by the control unit.

The integrated high temperature digester, with fast cooling feature, delivers precise and quick results. Even complex phosphate or nitrogen compounds are completely digested and processed in four to seven minutes; analogous to EN 1189 (phosphate) or ISO 11905-1 (nitrogen).

High quality, easy to handle GANICHEM reagents can be used for up to 100 determinations. Used reagents are safely and conveniently disposed of by Hach.

- Results in minutes
- Integrated fast digestion
- Ready-to-use reagents
- Portable control unit
- Simple handling

Part Number	Parameter	GBP Price
LPV401.99.00000	PO4-P	Contact Us