



Presence-Absence Media

Hach offers a variety of PA media to meet your testing needs. The modified Colitag method is the fastest P/A test on the market. PathoScreen products quickly screen for waterborne pathogens associated with fecal contamination. P/A products with MUG, compatible with paddle testers and BART products, are also available.

Part Number	Product description	GBP Price
2401650	Presence Absence Test with MUG, pk/50	Contact Us
2401612	E.coli and total coliforms (P/A) test, pk/12	Contact Us
2323250	Presence/absence test, pk/50	Contact Us
2323212	Presence/absence test, pk/12	Contact Us
2494925	Ampules, presence-absence broth, pk/25	Contact Us
2495525	E.coli and total coliforms, pk/25	Contact Us
2610696	PathoScreen Medium Presence/Absence Powder Pillows, pk/50	Contact Us
2610796	Pathoscreen PP 20 mL, pk/100	Contact Us
2401612-120	E.coli and total coliforms (P/A) test, 10x pk/12	Contact Us